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LABORATORY ROLLER KILN ER/P series

This type of kiln is made up of a sturdy steel structure repainted at 180 ° C with scratch-resistant epoxy paints.

Inside it encloses the thermal insulation consisting of low-density refractory bricks and pre-formed ceramic fiber plates highly resistant to temperature and thermal shocks.

The oven includes, in different quantities according to the various models, areas with different characteristics, namely:

- smoke entry and evacuation area
- preheating area
- firing zone
- cooling zone

Advancement system of the material:

The advancement of the tiles in the kiln takes place through ceramic rollers (with a diameter of 20 mm and a pitch of 38.1 mm), which are rotated by means of a pinion-chain system with programmable speed.

The towing system is developed to make it possible to disassemble the rollers in a sufficiently fast time.

Heating system

Depending on the model, the heating system consists of both spiral-wound Kanthal-type wire resistors and silicon carbide resistors, allowing these ovens to operate up to a maximum temperature of 1330 ° C.

They are inserted from the side of the oven allowing them to be easily replaced by removing only the protective casing; an operation that does not require the intervention of a specialized technician.

In the cooking zone, the upper and lower parts are managed independently. The cooling zones are managed automatically in all models, with the exception of the ER10 and ER15 models.

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
CONTROL PANEL

(12" TOUCH SCREEN terminal mounted on articulated support)

- 12" resistive TOUCH display, multilingual
- Storage of up to 30 recipes
- 3 types of recipe planning: manual / weekly / calendared
- Creation of a graph of recorded temperatures
- Privileged access to basic settings and advanced oven management settings
- Display bar and system error report
- ETHERNET port for remote control in teleassistance

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CONTROL PANEL

Contains the power section and all the equipment needed to control the furnace such as:

- main line switch
- motor overload cutouts
- fuses
- control boards
- electrical circuits for starting and controlling the various furnace devices
- static units

Kilns model	Max tile format [cm]	Firing cycle [minutes]		PROGRAMMABLES zones			NON-PROGRAMMABLES zones
		from	to	PREHEATING	FIRING	COOLING	COOLING
ER-15/P	15x15	11	110	1	1 (*)	/	1 (**)
ER-20/P	15x15	15	130	1	1 (*)	1	1 (***)
ER-25/P	20x20	20	174	1	1 (*)	1	1 (**) + 1 (***)
ER-30/P	25x25	13	125	2	1 (*)	2	1 (***)
ER-45/P	30x30	17	174	2	1 (*)	2	2 (**)

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
* : the temperature above and below the rollers can be individually programmed

** : the temperature is displayed on the control panel

*** : the temperature can be detected using an external portable thermocouple

The maximum size of the tile that can be fired is always related to the type of material and the duration of the firing cycle

TECHNICAL FEATURES

Mod.	Temp. max [°C]	Internal dimensions [mm]			External dimensions [mm]			Power kW	V + N	Weight [kG]
		[l]	[p]	[h]	[L]	[P]	[H]			
ER-15/P	1330	1500	250	40	2300	1100	1800	13	400	520
ER-20/P		2200			3400	1100	1800	16		680
ER-25/P		2630	350		3900	1200	1800	21		911
ER-30/P		3550	400		4750	1220	1800	36		1255
ER-45/P		4530			6160	1200	1800	40		1460