PRESSA DA LABORATORIO MIGNON S - MIGNON SS

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LABORATORY PRESS MIGNON S - MIGNON SS

Built in a compact and elegant structure of steel painted with epoxy paints cooked at 180 C°, they are fully hydraulic and are suitable for the production of samples in the laboratory.


The pressing force is exerted through a hydraulic cylinder placed on the upper part of the press and on which the appropriate punch is installed. The mold, placed on the work surface, is manually filled with the powder to be pressed and after pressing the sample obtained is extracted by a piston placed in the lower part.

The pressing force is manually adjustable and is displayed by means of a pressure gauge placed in the upper front part of the machine.

CONTROL PANEL

The front panel, to operate with the machine, presents:

- indicator light (shows if the machine is on or off)
- main switch
- two-handed control system
- extraction button

Supplied with:


- screed
- allen keys for mold assembly / disassembly

Not included:

- mold, supplied at the customer's choice

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Molde available on request:
 special molde Ø 40 mm.
 with sintered Zirconium Ox punch

NOTE

The moulds currently supplied are suitable for
 use with a maximum pressure on the shape of
 1700 Kg/cm²

IMPORTANT

For proper press operation, the pressure range of the circuit at
 which you can operate must be within between:

- 30 bar and 120 bar (MIGNON-S)
- 30 bar and 210 bar (MIGNON-SS)

Calculation example:

Piston diameter = 150 mm
 Piston area = 176,63 cm²
 Shape area = 50 cm²
 Desired pressure on the shape = 400 Kg/cm²
 Circuit pressure: $(400 \times 50) / 176,63 = 113,23$ bar
 (pressure value to be set)