PRESSA DA LABORATORIO MIGNON MAX68E PLUS

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LABORATORY PRESS MIGNON MAX68E PLUS

Built in a compact and elegant structure of steel painted with epoxy paints cooked at 180 ° C, it operates entirely hydraulically and is completely automatic except for the loading of the material to be pressed.

It is a suitable tool for the production of samples in the laboratory.

The pressing force is exerted through a hydraulic cylinder placed on the upper part of the press and on which the appropriate punch is installed.

The mold, placed on the work surface, is manually filled with the powder to be pressed and after pressing, the sample obtained is extracted by means of a piston placed in the lower part.

TECHNICAL SPECIFICATIONS

- Hydraulic operated pressing cylinder
- Hydraulically operated extractor cylinder
- Hydraulic power unit with automatic adjustment of the working pressure
- Electronic control unit that allow the press to work automatically
- Lateral aspirations to keep the work area clean
- Front lexan protection with safety microswitch
- Double hydraulic circuit

With this version of the press, it is possible to load a soft maximum of approximately 47 mm, therefore suitable for obtaining thick samples having a final thickness, when fired, of 20 mm.

The thickness of the soft can be programmed with a progressive value of 1 mm at a time. This means that with a single laboratory press we can have a soft height varying between 1 and 47 mm.

Therefore, this instrument is suitable for making pressed samples ranging from a few millimeters thick up to a maximum of 20 mm thick

ELECTRONIC CONTROL UNIT

- Visualization of the pressure on the display
- Programming of the working pressure
- Programming up to three thicknesses of soft 1
- Programming of the residence time at the maximum pressure set
- Programming of deaeration
- Start of the automatic pressing cycle

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
TOUCH SCREEN ELECTRONIC CONTROL UNIT

- Visualization of the pressure on the display
- Programming up to three thicknesses of soft 1
- Manual programming of the working pressure
- Programming of the residence time at the maximum pressure set
- Programming of deaeration
- Memorization of the mold dimensions
- Programming the Kg/cm2 desired in the pressing2
- display of the Kg/cm2 of force and relative working pressure, according to the format
- Start of the automatic pressing cycle
- predisposition with Ethernet and USB port for industry 4.0

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Mould available on request:

special molde Ø 40 mm. with sintered Zirconium Ox punch

NOTE

The moulds currently supplied are suitable for use with a maximum pressure on the shape of 1700 Kg/cm²

IMPORTANT

For proper press operation, the pressure range of the circuit at which you can operate must be within between 30 bars and 210 bars

Calculation example:

Piston diameter = 150 mm
Piston area= 176,63 cm²
Shape area = 50 cm²
Desired pressure on the shape = 400 Kg/cm²
Circuit pressure: $(400 \times 50)/176,63 = 113,23$ bar
(pressure value to be set)