PRESSA DALABORATORIO MIGNON MAX 40 E

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LABORATORY PRESS MIGNON MAX 40 E

Built in a compact and elegant steel structure coated with epoxy paint baked at 180°C.

It is entirely hydraulically operated and fully automatic except for loading the material to be pressed.

It is a tool designed for intermittent use, suitable for producing laboratory samples.

The pressing force is applied via a hydraulic cylinder located on the upper part of the press, to which a special pad is mounted.

The mold, placed on the work surface, is manually filled with the powder to be pressed. After pressing, the resulting sample is extracted by means of a piston located at the bottom.

MAXIMUM SOFT THICKNESS: 27 mm

TECHNICAL SPECIFICATIONS

- Hydraulically operated pressing cylinder (40 tons)
- Hydraulically operated extractor cylinder
- Hydraulic power unit with automatic adjustment of working pressure using a proportional valve (30-350 bar)
- Electronic control unit that allows the press to operate automatically
- Lexan front guard with safety microswitch
- 20 L hydraulic oil tank
- Bimetallic fixed-point safety thermostat 70°C
- Pressure transducer 0-400 bar
- Dual hydraulic circuit for independent calibration of specimen extraction force (factory setting)

TOUCH SCREEN ELECTRONIC CONTROL UNIT

- Pressure display on the display
- Programming of up to three soft material thicknesses (1)
- Manual programming of working pressure
- Programming of the dwell time at the maximum set pressure
- Prinogramming of deaeration
- Storing of mold dimensions
- Programming of the desired kg/cm2 pressing force (2)
- Display of the kg/cm2 force and the corresponding working pressure, depending on the format
- Automatic pressing cycle start
- Ready for Ethernet and USB ports for Industry 4.0


- (1) possibility of being able to press three types of powder simultaneously and therefore electronically regulate the thickness of the powders in the mould
- (2) the control unit will automatically adjust the working pressure

INCLUDED:

- leveler
- Allen keys for mold assembly/disassembly

MOULD NOT INCLUDED:

- supplied in a choice of standard or customized sizes with press-compatible features

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HYDRAULIC TANK:

- Oil level indicator/sight glass
- Oil filler cap
- Drain plug for emptying used oil/draining the tank

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Mod.	Largh.	Prof.	Alt.	Motore	V + N	Hz	Potenza	Pistone	Bar	Peso
	mm	mm	mm	kW			[ton]	[diam]	[max]	[kg]
MIGNON MAX 40 E	440	660	930	1,5	400	50/60	40	140	270	230

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On this model of laboratory press the moulds are interchangeable with each other and can be supplied in the following standard size:

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STANDARD SIZE MOLDS			
Size of the shape	case thickness	MIN on shape	MAX on shape
mm	mm	kg/cm ²	
Ø 40	30	368	1700
Ø 50		235	1700
30x80		192	1700
40x60		192	1700
50x50		185	1662
50x100		92	831
55x110		76	687
100x100		46	415
Other shapes, on request, according to compatibility			

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Mould available on request:

special molde Ø 40 mm. with sintered Zirconium Ox punch

NOTE

The moulds currently supplied are suitable for use with a maximum pressure on the shape of 1700 Kg/cm²

IMPORTANT

For proper press operation, the pressure range of the circuit at which you can operate must be within between 30 bars and 270 bars

Calculation example:

Piston diameter = 150 mm
Piston area = 176,63 cm²
Shape area = 50 cm²
Desired pressure on the shape = 400 Kg/cm²
Circuit pressure: $(400 \times 50) / 176,63 = 113,23$ bar
(pressure value to be set)